INSTRUCTIONS

tENS Vagal Boost

CONTENTS

What's in here.

01	Welcome	3
02	Meet Boost	4
03	In the box	5
04	Connecting the clip	6
05	Will it work with your phone?	7
06	If your phone won't power it	8
07	Before your first session	9
08	Your first fifteen minutes	10
09	The step everyone misses	11
10	The six modes	12
11	Finding your comfort level	13
12	Build the habit	14
13	If something feels off	15
14	Looking after Boost	17
15	Safety	18
16	Specs & warranty	19
17	Like Boost? Buddy & Oh!	20
—	About Vagal & getting help	21
—	Legal	22
—	Safety information (translated)	24

No battery. No charging. No app.

Boost is the simplest way to try vagus nerve stimulation. It plugs into your phone, borrows a trickle of power, and does the rest itself.

Your ear carries a branch of the vagus nerve — the nerve that tells your body *you're safe, you can relax now*. A gentle electrical pulse through a small clip on your tragus wakes it up. Fifteen minutes a day is enough.

Because Boost borrows its power, there is nothing inside to charge, pair, or update. Ever. It also means one thing really matters: *what you plug it into* — that's page 7, worth a minute before your first session.

THREE THINGS THAT TRIP PEOPLE UP

Hold M for two seconds to switch on — a quick tap does nothing.

Press + to at least level 1 or Boost turns itself off after thirty seconds.

Changing mode resets the level to zero — press + again every time.

Three buttons. That's the lot.

Display

Shows the mode letter (A–F) and, once set, the intensity level. Dark until you hold M — or if no power is arriving.

M button

Press and hold for about two seconds to switch on or off — a quick tap will not do it. Once on, short presses cycle through the six modes.

+ and –

Intensity, 1–9. Zeroes at power-on and each mode change — press + or nothing runs. Page 11.

Ear clip

Conductive-silicone electrode; the lead's pins push into it — page 6.

Cable

To a phone, laptop, or power bank.

03 • IN THE BOX

What you just unpacked.

Boost

The controller. Powered by whatever you plug it into.

Ear clip

Conductive-silicone auricular electrode.

Cable

Connects Boost to a USB-C port.

Lightning adapter

For iPhone 14 and earlier.

Comfort paste

A small pot. A thin layer on the pads softens the sensation.


Carry case

Everything lives here between sessions.

04 • CONNECTING THE CLIP

Pins into silicone.

Boost's lead ends in two small 2 mm pins. They push into the black silicone parts of the ear clip — and that's the whole connection.

Line up

One pin each side, into the black silicone sockets on the clip.

Push home

A firm, snug push is right — they're meant to hold.

To separate

Pull the pins, never the wires.

The black silicone *is* the electrode — its higher impedance is part of what keeps tENS gentle and safe.

05 • BEFORE YOU PLUG IN

Will it work with your phone?

Boost has no battery — your phone has to be willing to send power out of its port. Most are. Samsung phones are not.

✓ **WORKS**

iPhone 15 and later

Plug straight into the USB-C port. No adapter needed.

✓ **WORKS**

iPhone 14 and earlier

Use the Lightning adapter from the box, then carry on as normal.

✗ **NOT COMPATIBLE**

Samsung phones

Samsung handsets don't supply power to the Boost — across the whole range. No setting fixes it. See page 8.

✓ **WORKS**

Other Android

Motorola confirmed; most others plug straight in. On older models, switch on any USB OTG or power-output setting.

BATTERY POWER ONLY

Run Boost from a phone, laptop, or USB power bank — never from a mains wall charger, even to test it.

You haven't done anything wrong.

A dark screen almost always means the port isn't sending power — not that your Boost is dead.

IF YOU HAVE A SAMSUNG

That's the answer. We've seen it right across the range — A-series, S-series, Ultra, Flip, and tablets — and no setting changes it. Your Boost isn't faulty; your phone simply won't feed it. Older advice about a “USB OTG” toggle no longer applies: on modern Samsungs the setting doesn't exist, and where it does, it makes no difference.

THE PLAN THAT WORKS ANYWAY

Boost runs happily from a *laptop USB-C port* or a small *USB power bank*. A pocket power bank keeps the whole thing properly portable — which is rather the point.

RATHER RETURN IT?

Completely fair. Email help@vagal.com and we'll sort it without any fuss — you're inside the 45-day money-back guarantee.

What to expect.

Three habits make day one smooth.

HOLD M TO SWITCH ON

Press the M button *and keep holding it* for about two seconds, until the display lights up. A quick tap does nothing.

PREPARE THE PADS

Before you clip on, put a thin layer of the comfort paste from the box on both electrode pads. Dry pads concentrate the current and make it bite; paste spreads it evenly, and grips too. No paste? A drop of water works, or TENS gel from any pharmacy.

PRESS + STRAIGHT AWAY

Boost wakes at level 0 — no stimulation. Press + to at least level 1, or it switches off after thirty seconds. Every mode change resets the level to zero, so press + again.

Seven small steps.

1 Prepare the pads

Comfort paste on each pad, or a drop of water.

2 Clip the tragus

The small flap at the front of your ear. Left is the usual choice; either works. Pads on skin.

3 Plug in

Phone, laptop, or power bank. iPhone 14 or earlier — via the Lightning adapter.

4 Hold M (not a tap)

Press and hold for about two seconds — a quick tap does nothing. The display lights up once Boost is on.

5 Pick a mode

Short presses of M cycle A–F. Start with Mode C.

6 Press + to level 1

The step everyone misses. Boost wakes at level 0. If you don't press + within thirty seconds, it switches itself off. And every mode change resets the level to zero — so press + again.

7 Fifteen minutes

Boost won't stop on its own — hold M to finish. That's a session.

Turned on and chosen a mode? Now press +.

You held M for two seconds and the display lit up — good. Now power-on *and every mode change* reset the level to zero, and after thirty seconds at zero Boost switches itself off.

It isn't broken. It's waiting — and it resets on purpose: modes feel different, so a level that's comfortable in one could be sharp in another. The fix is always the same button: press +.

What you'll see

Press + and the level number appears next to the mode letter. Stimulation starts within a second or two.

Level 0

No stimulation, and switch-off thirty seconds later. Level 1 is the floor of the 1–9 scale — for most people it's close to the ceiling too.

DISPLAY GOES DARK MID-SESSION?

Changed modes and didn't press +? That's this. Level was set and it still cut out? See page 16.

10 • THE SIX MODES

A to F, decoded.

	PROTOCOL	IN BRIEF
A	Cardiac Coherence 38 Hz continuous	One steady rhythm, start to finish.
B	Multi-Frequency Cycle 22 → 10 → 5 Hz	Three rhythms in one session.
C	Low-Frequency 5 & 3 Hz alternating	The gentlest feel. Start here.
D	Autonomic Balance 38 & 10 Hz switching	High and mid, taking turns.
E	Dual High-Frequency 38 & 25 Hz switching	Two highs, taking turns.
F	Full Spectrum 78 → 38 → 25 Hz	The whole range in one session.

CHANGING RHYTHM IS NORMAL

Every mode except A changes frequency mid-session, with a brief pause at each switch — that's the protocol. And switching modes yourself resets the level to zero: press + again.

11 • FINDING YOUR COMFORT LEVEL

Comfortable beats intense.

Your skin doesn't need to hurt for the nerve to respond. If a session stings, one of these will fix it.

- | | |
|--------------------|--|
| 1. Contact | Comfort paste on the pads, every time. Water if you have none. |
| 2. Position | Clip centred on the tragus, both pads flat on skin. Off-centre contact concentrates the current and feels sharp. |
| 3. Level | Down a notch. Level 1–2 is where most people live. |
| 4. Mode | Not in Mode C? Switch to it — it's the gentlest. Then press + again: switching resets the level to zero. |

A TINGLE BEYOND THE EAR?

A mild echo along the jaw or neck can be normal at first — the nerve runs through that whole area. If anything is sharp, painful, or getting worse, stop and email us.

12 • BUILD THE HABIT


Fifteen minutes. Most days.

Some people feel calmer in their first session. The bigger differences build with steady daily use over four to six weeks — so give it that long before you judge it.

Anchor it to something you already do: the kettle, the evening news, the last emails of the day. Tick a box below each time. Missed a day? It doesn't reset anything — just start again tomorrow.

DAY

1	2	3	4	5	6	7
8	9	10	11	12	13	14

It won't light up.

First: **hold** M for two full seconds — a quick tap does nothing. Still dark? Then no power is arriving. It doesn't mean the device is dead, and it never needs charging.

SAMSUNG? THAT'S WHY

Samsung phones don't power the Boost. Skip straight to the test below.

THE TEST THAT SETTLES IT

Plug Boost into a **USB power bank** or a **laptop**, ear clip detached. Hold M for two seconds and watch the screen.

It lights up

Boost is fine — your phone is declining to send power. Use the power bank or laptop for sessions.

Still dark

That's a genuine fault. Email help@vagal.com and we'll replace it.

TWO THINGS TO RULE OUT FIRST

Keep your phone above 20% — many phones cut accessory power below that. And push the connector fully home; lint blocks more ports than faults do.

13 • IF SOMETHING FEELS OFF

Quick fixes.

NOTHING HAPPENS WHEN I PRESS M

Hold it — two full seconds. A quick tap does nothing. Still dark after holding? See page 15.

IT SWITCHES ITSELF OFF IN THIRTY SECONDS

Intensity is still at zero. Press + to at least level 1 — at power-on *and* after every mode change, because changing modes resets the level to zero. Page 11.

IT CUTS OUT MID-SESSION

Pressed M? That changed the mode and zeroed the level — thirty seconds later, off. Press + after every mode change. Didn't touch M? Try another battery source, keep your phone above 20%, and email us.

IT STINGS, OR I FELT A ZAP

Paste on the pads, centre the clip, drop a level, try Mode C. Still sharp at level 1? Stop and email us.

THE CLIP SLIPS OFF

Paste grips better than water. Using water? Less of it — damp, not wet.

I FEEL NOTHING / THE APP WON'T CONTROL IT

Feeling little at low levels is normal — it works either way. And there is no app: no Bluetooth, buttons only.

Almost nothing to do.

NOTHING TO CHARGE

There's no battery, so there's no charging — ever. If anyone tells you to charge it overnight, they're thinking of a different device.

AFTER EVERY SESSION

Wipe the clip electrodes with warm water, dry with a soft cloth, and put everything back in the case. Paste or gel left to dry is what wears pads out.

KEEP THE CONTROLLER DRY

Water on the electrodes: yes, a drop. Water on the controller: no. No taps, baths, or showers.

TRAVEL

No battery on board means nothing for airport security to mind. Boost goes anywhere your phone does — case in pocket, done.

SPARES

Replacement ear clips are at vag.al/shop if yours ever tarnishes or loses its spring.

Please read this page.

DO NOT USE BOOST IF YOU

Have a pacemaker or implanted electrical device. Are pregnant. Have a diagnosed heart rhythm problem. Experience persistent skin irritation at the stimulation site.

DO NOT USE BOOST WHILE

Driving or operating machinery. Near oxygen or flammable liquids. Showering, bathing, or swimming.

STOP AND SEE YOUR DOCTOR IF

You feel faint or unwell during or after a session. You're taking medication and unsure whether it interacts with vagus nerve stimulation.

A WELLNESS DEVICE, NOT A MEDICAL ONE

Boost is not intended to diagnose, treat, cure, or prevent any disease. For adults 18+ only.

The details.

STIMULATION

Method	Transcutaneous auricular VNS (tragus)
Modes	6 pre-set frequency programs
Frequency	3 Hz to 78 Hz
Waveform	Uniphasic square wave, 275 μ s pulse width
Output	Adjustable, low-intensity auricular — 9 levels
Electrodes	Conductive silicone; no metal at the skin

DEVICE

Power	USB-C, from a battery source — phone, laptop, or power bank. No internal battery.
Connectivity	Wired only. No Bluetooth, no app required.
Adapters	USB-C and Lightning included
Certification	CE, FCC

WARRANTY & SUPPORT

Twelve months against defects. Forty-five day money-back guarantee. help@vagal.com

17 • LIKE BOOST?

Try Buddy or Oh! next.

Boost has two bigger siblings, and everything you've learned here — placement, moisture, comfort — carries straight over.

Buddy carries its own battery and a joystick, and loads fresh session Tones from the app. Same clip, same ritual, nothing to plug into.

Oh! goes further still: fully wireless and in-ear, with a charging stand — tVNS that disappears into your day.

Curious? They're side by side at vag.al/compare — no hurry, Boost isn't going anywhere.

ABOUT VAGAL & GETTING HELP

A small team, actually here.

We're a small London-based team who've been building at-home tVNS since 2019 — first as Vagus.net, now as Vagal.com. Every device is built, tested, and supported by the same people. We read every email.

EMAIL US ANYTIME

help@vagal.com — we aim to reply within two business days. Check your spam folder if you haven't heard back.

This book, online

Scan for the always-current version of these instructions, or visit vagal.al/manual/boost

Ear placement guide: vagal.al/ear

The science

vagal.al/science

FAQs

vagal.al/faqs

LEGAL

The fine print.

TYPE BF

CE

FCC

PART 15

Type BF: body-floating protection against electric shock. CE compliant. FCC Part 15 certified.

LIMITED WARRANTY

Vagal Ltd warrants this Product against defects in materials and workmanship under normal use for twelve (12) months from the date of original purchase. This warranty does not cover damage caused by misuse, accident, abuse, external causes, use other than as described in this manual, finishes, normal wear and tear, tampering, or service by unauthorised agents. We do not warrant that the operation of the Product will be uninterrupted or error-free.

LIMITATION OF REMEDIES

Under this Limited Warranty, our liability and your exclusive remedy will be limited to replacement or repair of the Product. A replacement Product or part assumes the remaining warranty of the original Product or ninety (90) days from replacement, whichever is longer. To obtain warranty service, contact help@vagal.com. Return the Product with an insured carrier and retain your receipt. We are not responsible for items damaged or lost in transit. Other than for the reason of hardware defects, return freight cost is the customer's responsibility.

LEGAL (CONTINUED)

LIMITATION OF DAMAGES

In no event will Vagal Ltd or any affiliated company be responsible for any special, incidental, or consequential damages resulting from the use of this Product, or based on any breach of warranty, breach of contract, negligence, tort, or any other legal theory. Such damages may include, without limitation: loss of savings or revenue; loss of profit; loss of use; claims of third parties including retailers; cost of any substitute equipment or services. Some states do not allow the exclusion or limitation of incidental or consequential damages, so the above limitation may not apply to you. This warranty gives you specific legal rights, and you may have others which vary from state to state or country to country.

FCC STATEMENT

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation. Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

INTENDED USE

This device is intended for general wellness use. It is not intended to diagnose, treat, cure, or prevent any disease. If you have a medical condition, consult your doctor before use. Power only from a battery source (phone, laptop, or USB power bank); do not connect to a mains power adapter. For adults 18+ only.

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English UK · IE · MT

Do not use Vagal Boost if you have a pacemaker or implanted electrical device, are pregnant, have a diagnosed heart rhythm problem, or experience persistent skin irritation at the stimulation site. Do not use Vagal Boost while driving or operating machinery, near oxygen or flammable liquids, showering, bathing, or swimming. Stop and see your doctor if you feel faint or unwell during or after a session, or you're taking medication and unsure whether it interacts with vagus nerve stimulation. **A wellness device, not a medical one.** Vagal Boost is not intended to diagnose, treat, cure, or prevent any disease. For adults 18+ only.

Plug in. Clip on. Hold M. Press +.

Those last two matter.

Hold M for two seconds to switch on. Press + to start.

Tonight: Mode C, level 1, fifteen minutes.

Then most days for six weeks. That's the whole plan.

Better starts here.

Fifteen minutes a day, powered by the phone already in your pocket. The simplest way to look after your vagus nerve — and the friendliest price in tVNS.

Vagally Better®

Vagal Ltd
124 City Road
London EC1V 2NX
United Kingdom

help@vagal.com
vagal.com
US +1 456 VAGALLY
tENS Vagal Boost · VBtENS1