— INSTRUCTIONS

Wireless tVNS

CONTENTS

What's in here.

01	Welcome	3
02	Meet Oh!	4
03	In the box	5
04	The stand	6
05	Pairing	7
06	Fitting	8
07	Your first session	9
08	During the session	10
09	Tones	11
10	EverSteady	12
11	Building the habit	13
12	If something feels off	14
13	Looking after Oh!	15
14	Safety	16
15	Specifications	17
—	About Vagal & help	18
—	Legal	19
—	Safety information (translated)	21

You have the small one.

Eighteen grams. No screen, no buttons, no wires. Everything it does, it does quietly.

Your vagus nerve runs from your brainstem to your gut. It carries the signal that tells your body it is safe, and can stop bracing. When that signal is strong you sleep better, digest better, and recover from stress faster.

A branch of that nerve surfaces in your ear. It is the only place the vagus nerve can be reached from outside the body without surgery. A gentle current, delivered through a clip on your tragus, reaches it.

That is the whole idea. This book covers the rest: charging it, pairing it, wearing it, and what to expect once you do.

This is not a drug and not a surgery. It is a daily practice. The effects accumulate. Judge it at week six, not on day one.


No screen. No buttons.

Everything is controlled from your phone. There is no power switch — Oh! is awake whenever it has charge. Two lights are the only thing it tells you directly; battery, progress and intensity are all in the app.

LIGHT	MEANING
● Green, steady	Session running
● Green, fast blink	Still checking contact — reseal the clip
● Blue, slow blink	Connected to your phone, idle
○ White, steady	Charging
● Red, steady	Recovering a flat battery
● Red, fast blink	Updating firmware — leave it be
○ Off	Resting: not connected, not charging

What you just unpacked.

Oh!

The stimulator. 18 g, rechargeable, no buttons.

Charging stand

Holds Oh! while it charges. Cable runs in its channel.

Magnetic cable

USB-C. Seats into the top of the stand.

Ear hook

Silicone. Optional — see page 8.

Comfort Care

Conductive paste. Softens the sensation.

Cleaning cloth

For the electrodes after each session.

Tin

Keeps it all together between sessions.

Assemble this first.

The cable threads into the stand and stays there. Once only.

1 Seat the head

Push the magnetic head into the socket at the top of the stand.

2 Press the cable home

Hold the head in place and work the cable into the channel, down to the base.

3 Connect

Plug the USB-C end into any powered outlet, then rest Oh! on the head.

CHARGING

Sixty minutes of charge, sixty minutes of stimulation.

The app is the interface.

Oh! has no controls of its own. Without a phone it is an object. Pair it before your first session.

1

Install

Scan the code below, or search *Vagal Tones*. Free on both stores.

2

Find the device

With Oh! charged, open the Device page and let the app search. Once only.

3

Start

Clip Oh! to your ear, then press play. Intensity, length and mode are all set here.

vag.al/app">

vag.al/app

Scan to install. The page sends you to the right store for your phone.

Out of Bluetooth range the session pauses, and resumes when you come back.

Paste, clip, hook.

WHERE IT GOES

The tragus — the small flap of cartilage in front of your ear canal. Clip it so both electrodes touch skin.

COMFORT CARE PASTE


A smidge on both electrodes. It lowers the resistance between metal and skin, softening the sensation. Optional, but most prefer it.

THE EAR HOOK

Steadies Oh! if you are moving about, or if your ears are small. Sit the device inside the hook — rotating it within the hook changes the effective size.

WIDER FIT

NARROWER FIT

Five steps. Then wait.

1 Charge it

Sixty minutes on the stand from empty. The light tells you when it is done.

2 Paste both electrodes

A smidge of Comfort Care on each. Wipe your tragus with water first.

3 Clip it on

Onto the tragus, electrodes in contact. Add the ear hook if it needs steadying.

4 Start the session

Press play in the app. With *Auto* switched on, Oh! senses that it is seated and starts by itself.

5 Bring it up slowly

Raise until you notice it, then come back one step.

AIM FOR

Just below uncomfortable — not the highest you can stand. Higher is not better, and there is no prize for enduring it.

What normal feels like.

A TINGLE, NOT A SHOCK

A light prickling, like an itch from the inside. It usually settles within a minute as local blood flow rises and the skin adapts.

FEELING ALMOST NOTHING IS FINE

Especially with paste. Skin sensation and nerve stimulation are different things — the current is reaching the nerve whether or not you can feel it on the surface.

NOTHING DRAMATIC HAPPENS TODAY

This is not caffeine. There is no rush and no jolt. The effects are cumulative.

ADJUST OR STOP WHENEVER

Intensity is live in the app. Pause and resume as you like. If it becomes painful, stop — pain is not the goal and not a sign it is working.

WHEN IT ENDS

The current ramps down and the session closes. Unclip, wipe the electrodes, and put Oh! back on the stand.

One device. Many rhythms.

A Tone is a stimulation pattern: a frequency, a pulse width, a rhythm of on and off.

The vagus nerve does not respond to every pattern the same way. A slow, wide pulse does something different from a fast, narrow one — which is why Oh! is not fixed to a single setting.

Balance is already on your device. It runs at 30 Hz, it is the foundation Tone for daily use, and it is where everyone should start. You never have to buy anything to use Oh!

The rest of the library covers sleep, recovery, digestion, mood and sensitive skin. Each is bought once and syncs from the app.

vag.al/tones>

vag.al/tones

The current library, what each Tone is for, and the parameters behind it.

Changing Tone every day is the commonest way to get nowhere. Hold one for a fortnight, then change one thing.

EverSteady™

Your skin is a living organ. Its resistance shifts with heat, sweat and hydration.

Most stimulators push a fixed *voltage*, so when your skin changes the current changes with it: what you set in minute one is not what you get in minute twenty.

EverSteady regulates *current* instead. Fifty thousand times a second it measures the resistance, solves for the voltage it demands, and delivers the pulse you asked for.

50,000

CHECKS PER SECOND

1 μ A

RESOLUTION

10–60 V

SET AUTOMATICALLY

AUTO

Oh! can sense when it is seated correctly and start the session on its own. Switch it on in the app.

Adapt, Pulse+ and Switch: vag.al/eversteady

Fourteen days. That's the deal.

Consistency does more than duration. Fifteen minutes daily beats an hour on Sundays.

M	T	W	T	F	S	S
1	2	3	4	5	6	7
8	9	10	11	12	13	14

PICK A TIME YOU WILL ACTUALLY REPEAT

Morning coffee, the commute, the wind-down before bed. The specific hour matters far less than it being the same one.

MISSING A DAY IS NOT A RESET

Miss one, carry on. The fourteen days are a marker, not a test you can fail.

Almost always: contact.

Nearly every complaint traces back to the electrodes not meeting the skin properly.

SHARP OR STINGING

Too much current through too small a contact patch. Turn it down, add paste, and reseal the clip so both electrodes sit flat.

NOTHING AT ALL, AT ANY INTENSITY

Check the clip is on the tragus, not the lobe or rim. A fast-blinking green light means Oh! is still looking for good contact: reseal it and add paste.


IT CUTS OUT

Your phone has probably drifted out of Bluetooth range, or the battery has run down. Both resume where they left off.

REDNESS AFTERWARDS

Common and usually brief. If it persists between sessions, lower the intensity, use more paste, and give it a day. If it still persists, stop.

Small rituals. Long life.

WIPE THE ELECTRODES AFTER EVERY SESSION

Use the cloth. Paste left to dry on the contacts is the biggest cause of weak stimulation later. Cleaned each time, the electrodes last years.

DO NOT RINSE IT

Oh! is IPX4 — it shrugs off sweat and splashes, so workouts and rain are fine. It is not waterproof. Never submerge it, and never wear it in the shower or a pool.

KEEP IT ON THE STAND

Charge is the resting state. Between sessions, park it on the stand or in the tin.

THE BATTERY AGES

Like every lithium cell it holds less over the years. It does not fail suddenly; sessions just get shorter.

Please read this page.

DO NOT USE OH! IF YOU

Have a pacemaker or any implanted electrical device. Are pregnant. Have a diagnosed heart rhythm problem. Experience persistent skin irritation at the stimulation site.

DO NOT USE OH! WHILE

Driving or operating machinery. Near oxygen or flammable liquids. Showering, bathing, or swimming.

STOP AND SEE YOUR DOCTOR IF

You feel faint or unwell during or after a session. You're taking medication and unsure whether it interacts with vagus nerve stimulation.

ADULTS ONLY

For use by adults aged 18 and over.

A WELLNESS DEVICE, NOT A MEDICAL ONE

Vagal Oh! is not intended to diagnose, treat, cure, or prevent any disease.

The details.

Skin impedance drifts with temperature, moisture and movement. Oh! measures it 50,000 times a second and corrects every pulse, so the current reaching the nerve is the one you set.

STIMULATION

Frequency	0.1 Hz to 100 Hz (± 0.1 Hz)
Amplitude	0.1 mA to 4.0 mA (± 0.1 mA)
Max voltage	60 V
Pulse width	10 μ s to 5 ms (± 1 μ s)
Waveform	Biphasic square wave
Control	Current-controlled (EverSteady)
Impedance	Real time, 1–100 k Ω

DEVICE

Size, weight	23.5 × 23.5 × 19.6 mm, 18 g
Battery	86 mAh lithium polymer, 60 minutes
Charging	Magnetic stand, USB-C. 60 minutes.
Connectivity	Bluetooth Low Energy
Water	IPX4, splash resistant
Logging	8-hour session logs
Standards	IEC 60601-1, ISO 13485. CE, FCC.

Who made this thing.

We started in 2019 as Vagus.net, making electrodes for university research. We became Vagal.com in 2026. Tens of thousands of devices to more than forty countries, every one built, tested and supported by the same small team in London.

tVNS is not magic. It is a daily practice that quietly rebuilds a healthy nervous system. We do not promise cures. We build good tools, we say what the science actually shows, and we trust you to use them well.

GETTING HELP

help@vagal.com, or vagal.com/help for the contact form. A small team, considered answers, usually within two working days.

ONLINE

vagal.com/manual/oh for these steps with photographs. vagal.com/app for the app, protocols and firmware. vagal.com/faqs for questions we get often.

LEGAL

The fine print.

TYPE BF

CE

FCC

PART 15

Type BF: body-floating protection against electric shock (IEC 60601-1). CE compliant. FCC Part 15 certified.

LIMITED WARRANTY

Vagal Ltd warrants this Product against defects in materials and workmanship under normal use for twenty-four (24) months from the date of original purchase. This warranty does not cover damage caused by misuse, accident, abuse, external causes, use other than as described in this manual, finishes, normal wear and tear, tampering, service by unauthorised agents, or loss or damage to the battery. We do not warrant that the operation of the Product will be uninterrupted or error-free.

LIMITATION OF REMEDIES

Under this Limited Warranty, our liability and your exclusive remedy will be limited to replacement or repair of the Product. A replacement Product or part assumes the remaining warranty of the original Product or ninety (90) days from replacement, whichever is longer. To obtain warranty service, contact help@vagal.com. Return the Product with an insured carrier and retain your receipt. We are not responsible for items damaged or lost in transit. Other than for the reason of hardware defects, return freight cost is the customer's responsibility.

LEGAL (CONTINUED)

LIMITATION OF DAMAGES

In no event will Vagal Ltd or any affiliated company be responsible for any special, incidental, or consequential damages resulting from the use of this Product, or based on any breach of warranty, breach of contract, negligence, tort, or any other legal theory. Such damages may include, without limitation: loss of savings or revenue; loss of profit; loss of use; claims of third parties including retailers; cost of any substitute equipment or services. Some states do not allow the exclusion or limitation of incidental or consequential damages, so the above limitation may not apply to you. This warranty gives you specific legal rights, and you may have others which vary from state to state or country to country.

FCC STATEMENT

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation. This equipment complies with FCC radiation exposure limits set forth for an uncontrolled environment. Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

INTENDED USE

This device is intended for general wellness use. It is not intended to diagnose, treat, cure, or prevent any disease. If you have a medical condition, consult your doctor before use. For adults 18+ only.

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English UK · IE · MT

Do not use Vagal Oh! if you have a pacemaker or implanted electrical device, are pregnant, have a diagnosed heart rhythm problem, or experience persistent skin irritation at the stimulation site. Do not use Vagal Oh! while driving or operating machinery, near oxygen or flammable liquids, showering, bathing, or swimming. Stop and see your doctor if you feel faint or unwell during or after a session, or you're taking medication and unsure whether it interacts with vagus nerve stimulation. **A wellness device, not a medical one.** Vagal Oh! is not intended to diagnose, treat, cure, or prevent any disease. For adults 18+ only.

— FORTY-FIVE DAYS

Better, or your money back.

Over six weeks of daily use. Full refund if it isn't for you. We can offer that because most people who try it keep it.

[vag.al/manual/oh">](https://vag.al/manual/oh)

vag.al/manual/oh

These pages online,
with photographs and
anything added since
printing.

Vagal.com

Vagal Ltd
124 City Road
London EC1V 2NX
United Kingdom

Vagally Better®

help@vagal.com
vag.al/manual/oh
+1 456 VAGALLY
v0.8